Qufora® IrriSedo Bed


Qufora® IrriSedo **Bed**

EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Mode d'emploi	12
IT	Istruzioni per l'uso	16
ES	Instrucciones de uso	20
NL	Gebruiksaanwijzing	24
DA	Brugervejledning	28
SV	Bruksanvisning	32
NO	Bruksanvisning	36
HE	הוראות שימוש	40

Explaining the system


Explaining the Qufora IrriSedo Bed components

Water bag and pump


A water bag fitted with a valve and a pump.

Collection bag


A collection bag that holds 3.0 litres.

Rectal catheter

A rectal catheter with an insertion guide (for Ø21 and Ø30 catheters only).


System assembled


Instructions for use

To use the system please follow the instructions:


1. Fill the water bag with luke warm tap water (36–38° C) (only fill the bag up with the amount of water you need for your irrigation). Place the bag no higher than 50 cm above the level of the mattress. The length of the tube from water bag to the valve is 50 cm.
2. Fit the collection bag to the catheter with the insertion guide in the catheter. Note that the size Ø12 mm catheter does not require a guide.
3. Connect the tube to the rectal catheter.
4. Open the valve to fill the tube and pump with water. Close the valve again. Lubricate the tip of the catheter guide with a water based lubricant.
5. Hold the ring of the insertion guide with the scale of the collection bag pulled up between the 500 and 1000 ml. marking.
6. Position the patient on their left side. Carefully insert the rectal catheter into the rectum with the insertion guide.


1


2


3

4

5

6

7. With the catheter in the rectum pull out the insertion guide and let it drop into the collection bag.
 8. Push the plug in the collection bag into the open end of the catheter.
 9. Support the rectal catheter to avoid it being pushed out of the rectum. Open the valve on the tube. Start by opening the valve half way. If the water bag is placed at the mattress level or lower, open the valve completely and pump water into the rectum.
 10. Allow approximately 300 – 500 ml.(less for children – as specified by the healthcare professional) of water to enter the rectum and close the valve. Use max. 1000 ml. per irrigation. If any discomfort is felt, close the valve immediately.
- Please note: Should water seep out along the catheter the cause may be:**
- Too much water has entered the rectum.
 - The muscle tone of the sphincter is too low to hold the water.
 - The diameter of the catheter used is too small (The catheter is available in the diameters Ø12, Ø21 and Ø30 mm).
11. When the water is instilled, pull the plug out of the rectal catheter, water and faecal matter will flow into the collection bag.
 12. The procedure can be repeated by repeating step 1 to 4.
 13. Carefully remove the rectal catheter. Pull the water bag tube off the catheter. Empty the contents of the collection bag and dispose of the catheter and collection bag by placing them into a disposal bag before putting it into clinical waste. If required, all external surfaces of the system can be cleaned with lukewarm water and mild soap. Dry with a lint free cloth. Wash your hands to complete the procedure.

Instructions for use

Qufora IrriSedo Bed

Using the Qufora IrriSedo Bed makes it possible to perform transanal/rectal irrigation. Please read both pages of this instructions for use before using the Qufora IrriSedo Bed.

Indications for using transanal/rectal irrigation

The Qufora IrriSedo Bed is intended for adults and children 3 years and above to help manage faecal incontinence and/or constipation conditions where bowel management with rectal irrigation can only be undertaken when laying on the bed.

The Qufora® IrriSedo bed system requires the assistance of a carer to perform the procedure.

Use as advised by your healthcare professional

For the best results Qufora IrriSedo Bed should be used daily or as advised by your healthcare professional.

Before starting transanal/rectal irrigation

Use only after assessment and training from a qualified healthcare professional. A digital rectal examination is mandatory prior to your first irrigation. You should be supervised by an experienced healthcare professional during your first irrigation.

Transanal/rectal irrigation should not be used under the following circumstances (absolute contraindications):

- Known anal or colorectal stenosis
- Colorectal cancer/pelvic malignancy pre surgical removal
- Acute inflammatory bowel disease
- Acute diverticulitis
- Within 3 months of anal or colorectal surgery
- Within 4 weeks of endoscopic polypectomy
- Ischaemic colitis

Since the list is not exhaustive a healthcare professional should always consider individual patient factors as well.

Please note:

Transanal/rectal irrigation is not recommended for:

- Pregnancy or planned pregnancy or breastfeeding
- Children under 3 years of age

Use transanal/rectal irrigation only after careful discussion with a relevant medical practitioner under the following circumstances (relative contraindications):

- Inflammatory bowel disease (e.g. Crohn's disease or ulcerative colitis)
- Active perianal sepsis (fistula or abscess, third or fourth degree haemorrhoids)
- Previous rectal or colonic surgery
- Diarrhoea of unknown aetiology
- Faecal impaction/rectal constipation. If the patient is heavily constipated (faecally impacted) an initial clean-out of the bowel is mandatory before starting up on Qufora IrriSedo Bed Transanal Irrigation procedure.
- Severe autonomic dysreflexia
- Severe diverticulosis or diverticular abscess
- Abdominal or pelvic irradiation
- Long term steroid therapy
- Anti-coagulant therapy
- Low blood sodium
- Previous severe pelvic surgery
- Colonic biopsy within the past 3 months
- Use of rectal medications for other diseases which may be diluted by irrigation
- Congestive cardiac failure

Adverse events and device malfunction

Seek medical attention if the following is experienced during or after transanal/rectal irrigation:

- Sustained or severe bleeding from the rectum
- Sustained or severe abdominal or back pain, with or without fever

Perforation of the bowel wall is extremely rare. However it is a serious and potentially fatal complication of transanal/rectal irrigation. Should this occur, immediate medical attention is required. The incidents should be reported to the competent authority in your country and/or to Qufora A/S at QA@qufora.com. Qufora A/S should also be notified in case of any other device malfunction.

Please observe

Transanal/rectal irrigation can be associated with passing sensations of bowel discomfort, nausea, shivering, fatigue, sweating, headache and light bleeding from the rectum.

Please contact your healthcare professional for further advice.

Recommended usage time

The water bag and pump can be used up to 1 month. Should discolouration develop within the recommended usage time, this item should be replaced.

The catheter and collection bag are single use items.

If the single use items of the product are reused, there is a risk of unwanted effects that can influence the function of the product and the health of the patient, the carers and third persons.

Please find further information about this product on www.qufora.com.


Distributor/UKRP:

Qufora Ltd,
Euro House,
Satellite Park, Macmerrey,
East Lothian,
EH33 1RW

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