Qufora® IrriSedo MiniGo Flex

EN	Instructions for Use	2
DE	Gebrauchsanweisung	5
FR	Mode d'emploi	8
IT	Istruzioni per l'uso	11
ES	Instrucciones de uso	14
NL	Gebruiksaanwijzing	17
DA	Brugervejledning	20
SV	Bruksanvisning	23
NO	Bruksanvisning	26
FI	Käyttöohjeet	29
AR	لإمعة سالا تامي لعة	32
HE	שומיש תוארה	35

www.qufora.com

Qufora A/S
Gydevang 28-30
3450 Allerød
Denmark


Tel: +45 88 300 592
info@qufora.com
www.qufora.com

Explaining the system

Explaining the Qufora IrriSedo MiniGo Flex components

Soft pump


A soft pump (180 ml) with a flip lid that connects to the cone.

The pump can be used 30 times (equal to 1 month usage time*)

*when used for irrigation every day


Cone flex

A rectal cone with integrated flexible tube. The cone flex has a hydrophilic coating that is activated with water.

The cone flex is single-use


Assembled

The lid has a one-way valve, so the pump can be squeezed more than once, without air or water being sucked back into the pump during use. For this reason, the pump will stay compressed until the lid is opened.

Filling Qufora IrriSedo MiniGo Flex with water


1. Open the lid and fill the pump up to the inner rim with luke warm tap water (36-38° C) by holding it under the tap with the lid open.
2. Once the pump is filled to the inner rim hold onto the top part of the pump and press down on the lid until you hear the click sound.
3. Open the cone flex pouch at the perforated line. Whilst the cone flex is still in the pouch press the end of the tube firmly to the small round opening on the lid to assemble your system.
4. To activate the hydrophilic coating on the cone flex, wet the cone with water up to its base for a few seconds. Gently squeeze the pump until water comes out of the cone tip to remove the air. It is now ready to use.


1

2


3


4

Using Qufora IrriSedo MiniGo Flex

5. Sit on or over the toilet. Carefully guide the cone flex up to its base into the rectum. Now the cone is in place, gently squeeze the pump to instil water into the rectum. The pump can be squeezed several times to instil the water.
6. Remove the cone from the rectum. Faecal matter and water will be expelled into the toilet.


5

6

Emptying and disassembling Qufora IrriSedo MiniGo Flex



7. Place the used cone flex into the pouch and separate the two parts. The pouch can be folded and sealed using the sticky area, then disposed of in household waste. Please note that the cone flex is single-use only and must not be used again.

7

Cleaning and storing after use

- Open the lid of the pump to ensure that any residual water is emptied.
- If required, all external surfaces of the system can be cleaned with luke warm water and mild soap. Dry with a lint-free cloth.
- For storage leave the lid open.

Note: Please use Qufora IrriSedo MiniGo Flex as advised by your healthcare professional.

Instructions for use

Indications for using transanal/rectal irrigation

The Qufora IrriSedo MiniGo Flex is intended for adults and children 3 years and above to help manage faecal incontinence and/or constipation conditions.

For example:

- Passive faecal incontinence
- Post defaecation seepage
- Rectocele
- Incomplete evacuation
- Evacuation difficulties
- Low anterior resection syndrome

Within these groups some individuals may require the support of a carer to perform the procedure.

For the best results Qufora IrriSedo MiniGo Flex should be used daily or as advised by your healthcare professional.

Before starting transanal/rectal irrigation

Use only after assessment and training from a qualified healthcare professional. A digital rectal examination is mandatory prior to your first irrigation. You should be supervised by an experienced healthcare professional during your first irrigation.

Transanal/rectal irrigation should not be used under the following circumstances (absolute contraindications):

- Known anal or colorectal stenosis
- Colorectal cancer/pelvic malignancy pre surgical removal
- Acute inflammatory bowel disease
- Acute diverticulitis
- Within 3 months of anal or colorectal surgery
- Within 4 weeks of endoscopic polypectomy
- Ischaemic colitis

Since the list is not exhaustive a healthcare professional should always consider individual patient factors as well.

Use transanal/rectal irrigation only after careful discussion with relevant medical practitioner under the following circumstances (relative contraindications):

- Inflammatory bowel disease (e.g. Crohn's disease or ulcerative colitis)
- Active perianal sepsis (fistula or abscess, third or fourth degree haemorrhoids)
- Previous rectal or colonic surgery
- Diarrhoea of unknown aetiology
- Faecal impaction/rectal constipation
- Severe autonomic dysreflexia
- Severe diverticulosis or diverticular abscess
- Abdominal or pelvic irradiation
- Long term steroid therapy
- Anti-coagulant therapy
- Low blood sodium
- Previous severe pelvic surgery
- Colonic biopsy within the past 3 months
- Use of rectal medications for other diseases which may be diluted by irrigation
- Congestive cardiac failure
- Pregnancy or planned pregnancy or breastfeeding
- Children under 3 years of age

Adverse events and device malfunction

Seek medical attention if the following is experienced during or after transanal/rectal irrigation:

- Sustained or severe bleeding from the rectum
- Sustained or severe abdominal or back pain, with or without fever

Perforation of the bowel wall is extremely rare. However it is a serious and potentially fatal complication of transanal/rectal irrigation. Should this occur, immediate medical attention is required. The incidents should be reported to the competent authority in your country and/or to Qufora A/S at QA@qufora.com. Qufora A/S should also be notified in case of any other device malfunction.

Please observe

Transanal/rectal irrigation can be associated with passing sensations of bowel discomfort, nausea, shivering, fatigue, sweating, headache and light bleeding from the rectum.

Please contact your healthcare professional for further advice.

Recommended usage time

The soft pump can be used 30 times (equal to 1 month usage time*).

Should discolouration develop within the recommended usage time, the soft pump should be replaced.

The cone flex is for single-use only

If the cone flex is reused, there is a risk of adverse effects to the product and the health of the patient, carers and third persons.

Please find further information about this product on www.qufora.com.

*when used for irrigation every day.

Distributor/UKRP:

Qufora Ltd,
Euro House,
Satellite Park, Macmerry,
East Lothian,
EH33 1RW

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Qufora[®] is a registered trademark of Qufora A/S.© 02/2025. All rights reserved Qufora A/S, 3450 Allerød, Denmark.