
Qufora® IrriSedo **Klick**

Qufora® IrriSedo Klick


EN	Instructions for Use	4
DE	Gebrauchsanweisung	10
FR	Mode d'emploi	16
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	28
DA	Brugervejledning	34
SV	Bruksanvisning	40
NO	Bruksanvisning	46
AR	لإمعة سالإ تامي لعة	52
HE	הוראות שימוש	58

Explaining the system


System components


Rectal catheter
with a hydrophilic coating and
a water inflated balloon
- the catheter is single-use


Click tube
- can be used 15 times

Pump and control unit
with fixation strap
- can be used 90 times


Water bag - 1500 ml
with a tube and a valve
- can be used 15 times


Hook
to hang the water bag on
- can only be used on a clean
and smooth surface

System when assembled

The control unit

The system works by turning the control unit step by step in a 1-2-3-4 clockwise position.

- Position:**
- ① Use pump to prepare the system. Fill up tubes, pump and catheter with water, this empties the systems of air.
 - ② When the catheter is in the rectum, use the pump to inflate the balloon and create a water tight seal.
 - ③ Use pump or gravity to instill water into the rectum.
 - ④ The balloon will deflate and the catheter can be removed from the rectum.
 - ⌚ Turning the pointer to a **dot** will pause the system and stop the flow of water.


Instructions for use

Preparing the system for use

1. Wash and dry your hands.
2. Before you start, make sure the valve on the water bag tube is closed - lever is parallel to the tube.
3. Open the lid on the water bag and fill the water bag with the amount of luke warm tap water advised by your healthcare professional. Please note you need an extra 300 ml for preparing your system. Close the lid on the water bag.
4. Connect the water bag tube to the pump.
5. Connect the click tube to the control unit - when you hear a click it is connected.


1


2

3


4

5


Filling the system with water

6. Open the grey catheter pouch and connect the catheter to the click tube while leaving the catheter in the pouch - when you hear a click it is connected (6a). Take the catheter out without touching the hydrophilic coating (6b).
7. Attach the catheter pouch to a surface somewhere near the toilet, using the sticky area.
8. Turn the control unit to position 1.
9. Hold the pump with the control unit pointing upwards (9a) and open the valve on the water bag tube (9b). Squeeze the pump until it is full of water so that no air is left and water comes out of the catheter (9c).
10. Pause the system by turning the control unit to the dot position between 1 and 2.


Using the system


11. Wet the catheter up to its base with tap water. It is now lubricated and ready to use.
12. Hang or place the water bag no higher than shoulder height relative to your sitting position. You may use the supplied hook to hang the water bag on the wall. Now you are able to instill water by using gravity or pump.
13. Sit on or over the toilet. If required, fasten the control unit to your thigh or wheelchair, using the fixation straps.
14. Carefully insert the catheter into the rectum, supporting the catheter at its base.

15. Turn the control unit clockwise to position 2. Squeeze the pump once or twice (maximum of four complete pumps) to inflate the balloon with water (15a) until the catheter holds in the bowel by itself (15b). Your healthcare professional will advise you on this.
16. Turn the control unit clockwise to position 3. Water can now be instilled into the rectum. You can either pump the water into the bowel or you can let the water flow into the bowel using gravity. If you need to stop the flow turn the control unit to the dot between 3 and 4. To resume irrigation, turn the control unit back to position 3.
17. The amount of water flowing into the bowel can be monitored by using the scale on the water bag. Your healthcare professional will guide you regarding the amount of water to use for your irrigation. A general guide is 500 ml (this will vary for children) however this may be more or less depending on your requirement. If you feel any discomfort, do not instil more water into the bowel. Stop the flow by turning the control unit to the dot between 3 and 4.


18. When you have instilled the required amount of water into the bowel, turn the control unit clockwise to position 4. The balloon will now deflate and the catheter can be carefully removed from the rectum.
19. Close the valve on the water bag tube (lever is parallel to the tube).
20. Water and faecal matter will now be expelled. If this does not happen immediately, try coughing or leaning forward.
21. Click off the used catheter into the grey pouch by pressing the purple button. Hold the tube upwards to prevent water spillage. The pouch can be folded and sealed using the sticky area, then disposed of in household waste.


Emptying and disassembling the system


22. Place the system in the sink for emptying, separating and cleaning.
23. Separate the parts by clicking and pulling the connectors apart.
24. Empty the water from the system. To empty the control unit and the pump turn the control unit to position 1 and let the water run out through the control unit (24a). To empty the water bag open the valve on the water bag and let the water run through the tube (24b).


22

23

24a

24b

Cleaning and storing after use

25. If required, all external surfaces of the system can be cleaned with luke warm water and mild soap. Dry with a lint-free cloth.
26. Hang the system to dry with the water bag lid and the valve on the water bag tube open and leave the control unit in position 1.
27. Wash and dry your hands.

Instructions for use

Qufora Irrisedo Klick

Using the Qufora Irrisedo Klick makes it possible to perform transanal/rectal irrigation. Please read both pages of this instructions for use before using the Qufora Irrisedo Klick.

Indications for using transanal/rectal irrigation

The Qufora Irrisedo Klick is intended for adults and children 3 years and above to help manage faecal incontinence and/or constipation conditions.

For example:

- Urge faecal incontinence/urgency
- Constipation (slow transit/idiopathic/opioid induced/IBS-C)
- Neurogenic (spinal cord injury, upper/lower motor neurone/MS/parkinson's/spina bifida/cauda equina)

Within these groups some individuals may require the support of a carer to perform the procedure.

Use as advised by your healthcare professional

For the best results Qufora Irrisedo Klick should be used regularly as advised by your healthcare professional.

Before starting transanal/rectal irrigation

Use only after assessment and training from a qualified healthcare professional. A digital rectal examination is mandatory prior to your first irrigation. You should be supervised by an experienced healthcare professional during your first irrigation.

Transanal/rectal irrigation should not be used under the following circumstances (absolute contraindications):

- Known anal or colorectal stenosis
- Colorectal cancer/pelvic malignancy pre surgical removal
- Acute inflammatory bowel disease
- Acute diverticulitis
- Within 3 months of anal or colorectal surgery
- Within 4 weeks of endoscopic polypectomy
- Ischaemic colitis

Since the list is not exhaustive a healthcare professional should always consider individual patient factors as well.

Please note:

Transanal/rectal irrigation is not recommended for:

- Pregnancy or planned pregnancy or breastfeeding
- Children under 3 years of age

Use transanal/rectal irrigation only after careful discussion with a relevant medical practitioner under the following circumstances (relative contraindications):

- Inflammatory bowel disease (e.g. Crohn's disease or ulcerative colitis)
- Active perianal sepsis (fistula or abscess, third or fourth degree haemorrhoids)
- Previous rectal or colonic surgery
- Diarrhoea of unknown aetiology
- Faecal impaction/rectal constipation. If the patient is heavily constipated (faecally impacted) an initial clean-out of the bowel is mandatory before starting up on Qufora Irrisedo Klick Transanal Irrigation procedure.
- Severe autonomic dysreflexia
- Severe diverticulosis or diverticular abscess
- Abdominal or pelvic irradiation
- Long term steroid therapy
- Anti-coagulant therapy
- Low blood sodium
- Previous severe pelvic surgery
- Colonic biopsy within the past 3 months
- Use of rectal medications for other diseases which may be diluted by irrigation
- Congestive cardiac failure

Adverse events and device malfunction

Seek medical attention if the following is experienced during or after transanal/rectal irrigation:

- Sustained or severe bleeding from the rectum
- Sustained or severe abdominal or back pain, with or without fever

Perforation of the bowel wall is extremely rare. However it is a serious and potentially fatal complication of transanal/rectal irrigation. Should this occur, immediate medical attention is required. The incidents should be reported to the competent authority in your country and/or to Qufora A/S at QA@qufora.com. Qufora A/S should also be notified in case of any other device malfunction.

Please observe

Transanal/rectal irrigation can be associated with passing sensations of bowel discomfort, nausea, shivering, fatigue, sweating, headache and light bleeding from the rectum. Please contact your healthcare professional for further advice.

Recommended usage time

The water bag and the click tube can be used up to 15 times. The pump and control unit can be used up to 90 times.

Should discolouration develop on any part of the system within the recommended usage time, the part should be replaced.

The Klick catheters are single-use only

If the single-use items of the product are reused, there is a risk of unwanted effects that can influence the function of the product and the health of the patient, the carers and third persons.


Please find further information about this product on www.qufora.com.

Distributor/UKRP:

Qufora Ltd,
Euro House,
Satellite Park, Macmerry,
East Lothian,
EH33 1RW

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